

# Vanderbilt Particle Physics


**Professors Paul Sheldon, Will Johns,  
and Med Webster**

- **Global computing infrastructure** to support Particle Physics and other research making use of very large datasets (CMS at the LHC generates 0.5 Petabytes per year)
- **CMS detector** at the CERN Large Hadron Collider, general purpose detector, *discovered* the source of mass, the **Higgs** particle

# Cosmic Rays, high energy particles coming from space

The particles that hit our atmosphere:

- 90% protons (hydrogen nuclei)
- 9% alpha particles (helium nuclei)
- remainder solar wind particles and their *debris*



# Earth's Surface

The particles that hit our atmosphere clobber the nuclei of the atoms and molecules:

- create a shower of x-rays, gamma rays, and light
- also make a lot of muons (like an electron but heavier), with mass 106 MeV (compare to an electron with mass 0.511 MeV).


# Energy and Quantity

**eV** is short for electron-Volt,  
and gives the energy of an  
object...or its mass!

$$E = mc^2$$

**Energy equals mass times  
speed of light squared** (just  
some constant, but a big one)

proton mass 0.924 GeV  
electron mass 0.511 MeV  
muon mass 106 MeV  
my mass  $4.3 \times 10^{37}$  eV



# Speed of Light

Fastest possible speed is the speed of light in vacuum.

**Defined as**  $299792458 \text{ m/s}$

$$3.0 \times 10^8 \text{ m/s}$$

$$30 \text{ cm/ns}$$

$$300 \text{ m}/\mu\text{s}$$

$$300 \mu\text{m}/\text{ps}$$

# SI Prefixes


**Table 5. SI prefixes**

Factor	Name	Symbol	Factor	Name	Symbol
$10^{24}$	yotta	Y	$10^{-1}$	deci	d
$10^{21}$	zetta	Z	$10^{-2}$	centi	c
$10^{18}$	exa	E	$10^{-3}$	milli	m
$10^{15}$	peta	P	$10^{-6}$	micro	$\mu$
$10^{12}$	tera	T	$10^{-9}$	nano	n
$10^9$	giga	G	$10^{-12}$	pico	p
$10^6$	mega	M	$10^{-15}$	femto	f
$10^3$	kilo	k	$10^{-18}$	atto	a
$10^2$	hecto	h	$10^{-21}$	zepto	z
$10^1$	deka	da	$10^{-24}$	yocto	y



# QuarkNet Detectors

Flash description at  
[http://www18.i2u2.org/elab/cosmic/flash/daq\\_only\\_standalone.html](http://www18.i2u2.org/elab/cosmic/flash/daq_only_standalone.html) .






# Pieces / Parts

**Scintillation Counter:** Cosmic Rays create light in the scintillators. The photomultiplier converts light to an electronic pulse.


**Data Acquisition (DAQ) Board:** When the signal is in a certain range, the DAQ executes onboard logic and processes the signal. The result may be a string of data sent to a local computer.



**Computer:** The DAQ sends data to a local computer for storage and later analysis. The hexadecimal text includes information on pulse size and GPS time for signals from each counter.

# Detectors Details 1/3

DAQ board is **Quarknet v2.5**. It has 4 coax connectors for the PMT, and ethernet for the GPS, USB to communicate with the computer, input power, and a power out to an electrical box that powers the 4 PMTs. PMT voltages are no. 1, 550V; no. 2, 735V; no. 3, 925V, and no. 4, 885V, approximately.



Coax from PMT

GPS Conn

Power IN and OUT  
to Power Cond. for  
the PMTs

USB to Computer

# Detectors Details 2/3


Scintillator "paddle" is 12x10 inches connected to a cylindrical photomultiplier tube (PMT).





# Detectors Details 3/3

On the computer a program written in LabView configures the DAQ card, receives the data, and writes to the computer files.



# Data Handling - Spreadsheet 1/4

You must take data in the **Excel tab**, presumably a file with a CSV extension. The file is actually "tab separated." In an editor...



time	r0	f0	r1	f1	r2	f2	r3	f3	
62325.050781		-1	-1	-1	-1	29	38	28	39
62355.062500		-1	-1	29	49	29	46	31	40
62550.562500		-1	-1	22	46	20	35	17	42
62616.429687		-1	-1	-1	-1	-1	-1	-1	-1
62738.121094		-1	-1	17	35	11	29	12	23

**time** - the time stamp in units of seconds

**r0, r1, r2, r3** - the time at the rising edge of the corresponding PMT pulse

**f0, f1, f2, f3** - the time at the falling edge of the corresponding PMT pulse

**rN** and **fN** time is in "ticks," **40/32 ns = 1 tick**. Older boards use 24/32 ns = 1 tick. The two boards use two different internal clocks, new at 25.0 and older at 41.666 MHz.

Google Drive folder, <http://drive.google.com>

# Data Handling - Spreadsheet 2/4

- Copy the file, xxxx.csv into your Google Docs/Drive area.
- Rename it xxxx.tsv (for tab separated values).
- Right Click and open the file the Google Sheets.
- Extend out the columns for our calculations: click on column “I” (the last one) and extend the column 1 right. You can use the down arrow that shows in the column. Add about 6 columns.
- Leave one blank column between the left table of numbers and the set we will calculate.
- In the next column over, “K”, put a title like “r3-r1” (that is the difference in the risetime for paddle no. 4 and no. 2).
- Under that column create a math formula for the difference, sometime like “=h2-d2” and then Copy and Paste for the whole column. Check a couple of the numbers by eye/hand. Make sure it is working.
- If there were no other delays in the signals to the computer, that column would give the speed of the particles through the paddles.

# Data Handling - Spreadsheet 3/4

- Start looking at the differences of the risetimes.
- The units on the rise and falltimes is  $40/32 \text{ ns} = 1.25 \text{ ns}$ .
- Remember the speed of light in useful units is  $30 \text{ cm/ns}$ .
- For the list of measurements of interest, do some simple statistics in the spreadsheet:
  - Count the number, “= count(k2:k30)”, find the average or mean, “=average(k2:k30)”, and the standard deviation or spread in the measurements (not the *error in the mean*), “=stdev(k2:k30)”.



# Data Handling - Spreadsheet 4/4

- Particle Physicists love to plot "histograms" or frequency plots. In Spreadsheets you use the "frequency" function to make a column of frequency data and then you plot it as a Scatter plot.
- It looks something like "`=frequency(k2:k30; n2:n21)`". The data is k2:k30 and the set of numbers to use as bins are n2:n21.

	I	J	K	L	M	N	O	P	Q	R	S	T
f3			r3-r2	r1-r0	r2-r0	r3-r1						
7	33		4	-4	4	12	14 max					
4	35		3	-3	2	8	-10 min					
							from min to max					
5	25		5	0	2	7		frequen	bins	frequen		
9	69		0	-3	5	8	-10	1	0	1		
2	35		6	-1	1	8	-8	0	1	0		
3	43		0	11	19	8	-6	0	2	0		
5	28		5	0	4	9	-4	0	3	0		
5	45		-3	-5	4	6	-2	0	4	0		
5	52		3	-8	0	11	0	0	5	0		
3	63		0	3	8	11	2	3	6	3		
0							4	19	7	19		
3							6	44	8	44		
3							8	73	9	73		
3							10	49	10	49		
9							12	23	11	23		
5							14	3	12	3		
7							16	0	13	0		
2								0	14	0		
3									15			
4												
2												
0												
9												
3												
3												
4												
1												

In top cell, =frequency(N2:N216)

peak at  
ns/tick times      9 ticks  
                        0.75 ns/tick  
                        6.75 ns

**In top cell, =frequency(N2:N216,O5:O18)**

peak at	9 ticks
ns/tick	0.75 ns/tick
times	6.75 ns

# Data Handling - Results

- Look for changes in time that are all about the same.
- Toss out any obviously bad data
- Remember average speed is "change in distance / change in time"

$$V_{avg} = \frac{\Delta D}{\Delta t}$$

# Websites:

Vanderbilt Particle Physics, <http://www.hep.vanderbilt.edu/~focuser/index.php>

Particle Adventure, <http://www.particleadventure.org>

Interactions (news), <http://www.interactions.org/cms>

Compact Muon Solenoid at the CERN LHC, <http://cms.cern.ch>

Cosmic Rays (wikipedia), [http://en.wikipedia.org/wiki/Cosmic\\_ray](http://en.wikipedia.org/wiki/Cosmic_ray)

QuarkNet, <http://quarknet.fnal.gov>

Cosmic Ray e-Lab, <http://www18.i2u2.org/elab/cosmic/home/>

Detector Description, [http://www18.i2u2.org/elab/cosmic/flash/daq\\_only\\_standalone.html](http://www18.i2u2.org/elab/cosmic/flash/daq_only_standalone.html)

New Manual for DAQ Card